

# Datasheet

## Excel Life™ Medical Range, Common Plate Horizontal Single Powerpoint 10A – Antimicrobial\*

with Power Available LED, marked “RCD PROTECTED” and Circuit ID option

CAT. NO. EC787NCAWE, EC787NCABU, EC787NCARD



EC787NCAWE



EC787NCABU



EC787NCARD


Cat. No.	EC787NCAWE, EC787NCABU, EC787NCARD
Colour	White(WE), Blue(BU) or Red(RD)
Description	Excel Life™ Medical Range Common Plate – Horizontal Single Powerpoint 10A – Antimicrobial with Power Available LED and marked “RCD PROTECTED” (Circuit ID option)
Nominal Input Voltage	230 - 240 V a.c. 50 Hz
Current	10 A
Protection against solid bodies	IK04
IP rating	IP30
Power Available LED	Amber
Plate size (W x H x D)	118 x 74 x 10 mm
Mounting centres	84 mm
Operating temperature	5° to 35°C
Glass/tile cut out	56 x 96 mm
Compliance	AS/NZS 3112
Product Marking	

# Excel Life™ Medical Range, Common Plate Horizontal Single Powerpoint 10A – Antimicrobial\*



with Power Available LED, marked "RCD PROTECTED" and Circuit ID option

CAT. NO. **EC787NCAWE, EC787NCABU, EC787NCARD**



## FEATURES

- Chemical Resistant and Antimicrobial
- Internal safety shutters for extra safety
- Power Available LED (Amber)
- More fixing options – 4 corners
- Retained terminal screws
- Available with Circuit ID
- Marked with "RCD protected"
- Colour White, Blue (UPS), Red (Essential)

Single product		
Cat. No.	Barcode	Dimension (W x H x D)
EC787NCAWE	9321001427841	85 x 117 x 50 mm
EC787NCABU	9321001427858	
EC787NCARD	9321001427865	

Inner box (QTY 10)		
Cat. No.	Barcode	Dimension (W x H x D)
EC787NCAWE	9321001428657	175 x 165 x 125 mm
EC787NCABU	9321001428664	
EC787NCARD	9321001428671	

Outer box (QTY 120)		
Cat. No.	Barcode	Dimension (W x H x D)
EC787NCAWE	9321001429463	608 x 308 x 340 mm
EC787NCABU	9321001429470	
EC787NCARD	9321001429487	

*Antimicrobial Specifications			
<b>Maintenance</b>	Surface cleaning with cloth: Hospital-grade water-based cleaners only. Do not use any solvent based cleaner such as acetone, tar-remover or trichloroethylene.		
<b>Antimicrobial Properties</b>	Challenge Organism	Polycarbonate alloy in fascia parts	Polycarbonate alloy in switch dolly
		Reduction Coefficient/Kill Rate	
<b>Note:</b> 5 Log = 99.9999% 4 Log = 99.99% 3 Log = 99.9% 2 Log = 99%	<b>Escherichia Coli</b>	2.55 Log / 99.72%	6.72 Log / 99.9999%
	<b>Staphylococcus aureus</b>	2.70 Log / 99.79%	5.22 Log / 99.999%
	Tests performed according to JIS Z 2801:2010 (Antimicrobial products – Test for antimicrobial activity and efficacy). Test method is equivalent to procedure detailed in ISO 22196:2011 (Measurement of antibacterial activity on plastics and other non-porous surfaces). Test results available upon request. Not to be used as an alternative or substitute for standard hygienic practices, procedures, or protocols. This device has a long-term antimicrobial effect. Lifespan is dependent upon cleaning procedures and the use of recommended cleaning products. Product manufactured with antimicrobial material based on Ag+ ions. This technology limits the proliferation of bacteria on our devices without increasing antibiotic resistance.		

Legrand Australia  
1300 369 777  
www.legrand.com.au  
ABN 31 000 102 661

Legrand New Zealand  
0800 476 009  
www.legrand.co.nz

Refer to your group buying office, Legrand Sales Representative or Legrand Sales Office for pricing.

Please note: specifications are subject to change without prior notice.



LE15183AA

08/2024